

# Curved Femoral Head Impactor

Designed by Amiee Zirpel

*Allows for in-line femoral head impaction during minimally invasive THR*

The curved offset handle allows the head impactor to be slid under the skin of a small incision, and helps provide hand-held stability and maneuverability within the wound, while the impaction platform is easily accessible outside the wound. The impaction disc is made of delrin, which helps prevent marring and scratching of components.

**PRODUCT NO:**

3644

Overall Length: 7.25"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

**INNOMED**

1.800.548.2362