


# DAA Canal Finder Rasp

*Designed to help begin preparation of the femoral canal prior to stem broaching – features a large handle with a striking plate end*



**PRODUCT NO:**

C1026

Overall Length: 13.125" (33,4 cm)

Handle Length: 6" (15,2 cm)

Maximum Shaft Diameter: 13 mm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**