

# Direct Anterior Approach Retractor Set

A General Use Set of Innomed Retractors for Direct Anterior Approach Total Hip Arthroplasty

Single Prong Acetabular Retractor – Standard  
#6570



Mueller-type Femoral Neck Elevator – Standard  
#3415



Cobra Retractor – Standard with Sharp Tip  
#6129



Deep Hohmann-style Retractor with Large Handle – Standard  
#C1009



Modified Hohmann Retractor – Narrow  
#4535



Cobra Retractor – Narrow  
#6120



Bent Hohmann Retractor – Narrow with Extra Long Handle  
#7110-01



Bone Hook – Large  
#5920



Rivero Extra Grip Femoral Head Remover with Zimmer Hall Quick-connect  
#3706



PRODUCT NO'S:
6500-01 [DAA Retractor Set] Set includes (2) #6120 and (1) of the other instruments shown below
<b>Set Includes / Available Individually:</b>
6570 [Single Prong Acetabular Retractor – Standard] Overall Length: 12.125" (30,8 cm) Blade Width: 22.3 mm
3415 [Mueller-type Femoral Neck Elevator – Standard] Overall Length: 13.5" (34,3 cm) Handle Length: 6.5" (16,5 cm) Blade Width at Widest: 25 mm
6129 [Cobra Retractor – Standard w/Sharp Tip] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade at Widest: 33 mm
C1009 [Deep Hohmann-style Retractors with Large Handle – Standard] Overall Length: 15.75" (40 cm) Handle Length: 6.125" (15,6 cm) Blade Depth: 4.5" (11,4 cm) Blade Width: 1.1" (2,8 cm) Prong Length: 1.5" (3,8 cm)
4535 [Modified Hohmann Retractor – Narrow] Overall Length: 10" (25,4 cm) Blade Width: 14 mm
6120 [Cobra Retractor – Narrow] <b>(2) included in set, (1) only with this product number</b> Overall Length: 11.75" (29,8 cm) Handle Length: 6.5" (16,5 cm) Blade Width: 19 mm
7110-01 [Bent Hohmann Retractor – Narrow with Extra Long Handle] Overall Length: 11.5" (29,2 cm) Handle Length: 10" (25,4 cm) Blade Width: 19 mm Depth from Bend: 4.25" (10,8 cm)
5920 [Bone Hook – Large] Curve Diameter: 50 mm Overall Length: 12.75" (32,4 cm) Handle Length: 4.75" (12,1 cm)
3706 [Rivero Extra Grip Femoral Head Remover with Zimmer Hall Quick-connect] Overall Length: 8.5" (21,6 cm)



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.