


Deep Hohmann-style Retractors with Large Handle

Designed for retraction around the femur and acetabulum



PRODUCT NO:

C1009 [Standard]

Overall Length: 15.75" (40 cm)

Handle Length: 6.125" (15,6 cm)

Blade Depth: 4.5" (11,4 cm)

Blade Width: 1.1" (28 mm)



PRODUCT NO:

C1010 [90°]

Overall Length: 14.25" (36,2 cm)

Handle Length: 6.125" (15,6 cm)

Blade Depth: 5.25" (13,3 cm)

Blade Width: 1.1" (28 mm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362