



# Delrin Insert Pliers

*Designed to grasp an implant for adjustment without marring the implant surface*



**PRODUCT NO'S:**

2025

*Overall Length: 8*

2025-03 [Replacement Insert]

*Includes top and bottom delrin jaws,  
two screws and a hex wrench*

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



**INNOMED** 

**1.800.548.2362**