


## Desai Curette Osteotome

Designed by Sarang Desai, DO

*Designed to remove bone and cartilage, helpful for preparing joint surfaces for fusion, allowing easy removal of osteophytes and cartilage without having to switch instruments*

The osteotome portion also can be used to "feather" the subchondral surface to expose bleeding bone. It is also useful in instances of obtaining autograft, as it can be used to create a bone window and then remove cancellous bone.



**PRODUCT NO:**  
5241  
Overall Length: 8.25" (21 cm)  
Handle Length: 4.25" (10,8 mm)  
Cup: 3.5 x 5 mm  
Tongue Width: 3,5 mm  
Tongue Length: 3.5 mm from edge of cup



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED** 

1.800.548.2362