

# Desai Surgical Funnel


Designed by Sarang Desai, DO

*Helps with control and placement of bone graft or antibiotic beads*

Made from surgical grade stainless steel (for sterilization purposes).



**PRODUCT NO:**  
8989  
Overall Length: 6.25" (15,9 cm)  
Handle Length: 3.25" (8,3 cm)  
Funnel Diameter at Top: 3" (7,6 cm)  
Funnel Flow-thru Diameter: 11 mm



Profile View

ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

# INNOMED

**1.800.548.2362**