

Disposable LED Light Source

Compact, disposable source of brilliant 6000K white light for use with lighted retractors



Single Use, Sterile, LED Lightsource

- Compact, disposable source of brilliant 6000K white light.
- For use with lighted retractors and rigid or flexible endoscopes.
- Connects directly to device via integral ACMI light port.
- Enables untethered operation for easier positioning.
- Eliminates cost of reprocessing a fiberoptic cable.
- Light output approximately 70 lumens.
- Simple pull-tab activation; automatic shutoff at 1.5 hrs. nominal.
- Packaged individually in EtO sterile pouch; 10 pouches per box.

Clinical Benefits:

- Eliminate chance of cross-contamination.
- Cost savings by not having to reprocess through decontamination.
- Save time and improve efficiency in the OR.

Essentials:

- Dimensions: 2.5" (6.4cm) long; 1.12" (2.8cm) diameter.
- Weight: 1.1 oz. (31g).
- FDA Class I device, **CE** marked.
- RoHS, REACH and WEEE compliant.



PRODUCT NO'S:
PACKAGE OF 1:
8010-01 [Disposable LED Light Source] Overall Length: 2.5" (6.4 cm) Diameter: 1" (2,54 cm)
PACKAGE OF 10:
8010-10 [Disposable LED Light Source]



ISO 13485:2016

Scan to
Launch Our
Website



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362