



# Dodson Extremity Skin Saver

Designed by Mark A. Dodson, MD

*Designed to help protect the patient's skin  
when removing a disposable tourniquet*

After the device is laid on the skin, the disposable tourniquet is then rolled over it, and stays in position until removed at the end of the case.

**PRODUCT NO:**

8628

Overall Length: 4.75" (12,1cm)

Width: 1.5" (3,8cm)

Lip: .5" (1,3cm)



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71



**1.800.548.2362**