

# Dodson Extremity Skin Saver

Designed by Mark A. Dodson, MD

*Designed to help protect the patient's skin  
when removing a disposable tourniquet*

After the device is laid on the skin, the disposable tourniquet is then rolled over it, and stays in position until removed at the end of the case.




**PRODUCT NO:**

8628

Overall Length: 4.75" (12,1 cm)

Width: 1.5" (3,8 cm)

Lip: .5" (1,3 cm)



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**