



Double Ended Grater Cleaning Tool

Designed by Brandon Thompson, CST/CFA

Designed for right or left handed use to easily remove bone fragments from acetabular graters

PRODUCT NO:

8007

Overall Length: 7" (17,8 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362