

Duke Classic Acetabular Retractor with Extra Grip Tip

Designed by Justin Duke, MD

*Designed to retract the femur during acetabular exposure
for either posterior or lateral approaches*

*Posterior approach usage demonstrated. May also be used
to retract the femur posteriorly during lateral exposure.*




PRODUCT NO:

7622
Overall Length: 13.5" (32,3 cm)
Handle Length: 9" (22,9 cm)
Depth from Bend: 4" (10,2 cm)
Prong Depth: 1.25" (3,2 cm)



*Posterior approach usage
demonstrated. May also be used
to retract the femur posteriorly
during lateral exposure.*



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362