



## Duke Anterior/Posterior Inferior Retractors with Extra Grip Tip

Designed by Justin Duke, MD

*An inferior acetabular retractor designed for total hip arthroplasty while prepping the acetabulum*

### PRODUCT NO'S:

7621-01 [Left]  
Overall Length: 10.25" (26 cm)  
Handle Length: 5" (12,7 cm)  
Depth: 7.5" (19,1 cm)  
Blade Length: 2.75" (7 cm)  
Blade Width: 1.125" (2,9 cm)  
Prong Length: 1" (2,5 cm)  
Prong Width: 6 mm

7621-02 [Right]  
Overall Length: 10.25" (26 cm)  
Handle Length: 5" (12,7 cm)  
Depth: 7.5" (19,1 cm)  
Blade Length: 2.75" (7 cm)  
Blade Width: 1.125" (2,9 cm)  
Prong Length: 1" (2,5 cm)  
Prong Width: 6 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED** 

1.800.548.2362