

Duke Classic Inferior Retractors with Extra Grip Tip - Left & Right


Designed by Justin Duke, MD

An inferior acetabular retractor designed for total hip arthroplasty while prepping the acetabulum

PRODUCT NO'S:

7621-01 [Left]
Overall Length: 10.25" (26 cm)
Handle Length: 5" (12.7 cm)
Depth: 7.5" (19.1 cm)
Blade Length: 2.75" (7 cm)
Blade Width: 1.125" (2.9 cm)
Prong Length: 1" (2.5 cm)
Prong Width: 6 mm

7621-02 [Right]
Overall Length: 10.25" (26 cm)
Handle Length: 5" (12.7 cm)
Depth: 7.5" (19.1 cm)
Blade Length: 2.75" (7 cm)
Blade Width: 1.125" (2.9 cm)
Prong Length: 1" (2.5 cm)
Prong Width: 6 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362