

Durham Offset Kolbel Shoulder Retractor Set

Designed by Alfred A. Durham, MD

Designed for retraction of the deltoid and under the short head of the biceps muscle to expose the shoulder, the longer offset blades are useful in patients with large muscles, and the shorter offset blades are useful in smaller elderly patients

Snap-in, freely pivoting smooth curved blades help to concentrate the forces on the center of the muscle bellies, allowing the retractor to remain centered and not get in the way of exposure.



Set comes with retractor handle (T1030-01) and 1 pair each of the Long Offset Blades (T1030-L) and the Short Offset Blades (T1030-S).

PRODUCT NO'S:
T1030 [Set] Set includes: (1) T1030-01, (2) T1030-L, (2) T1030-S
Also available individually:
T1030-01 [Retractor Handle] Overall Length: 7" (17,8 cm)
T1030-L [Long Offset Blade] (2) included in set, (1) only with this product number Offset Length: 35 mm Blade Dimensions: 36 x 36 mm
T1030-S [Short Offset Blade] (2) included in set, (1) only with this product number Offset Length: 10 mm Blade Dimensions: 36 x 36 mm



36 x 36 mm
10 mm Offset

36 x 36 mm
35 mm Offset

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362