

Eickmann Knee Revision Set

Designed by Thomas Eickmann, MD

Used for total knee revision



PRODUCT NO'S:

5470-00 [Complete Set]

Individual Instruments:

5470-08 [8 mm Chisel]

Osteotome Width: 8 mm
Blade Length: 2.375" (6 cm)
Overall Length: 7.375" (18,7 cm)

5470-11 [11 mm Chisel]

Osteotome Width: 11 mm
Blade Length: 2.375" (6 cm)
Overall Length: 7.375" (18,7 cm)

5470-20 [20 mm Chisel]

Osteotome Width: 20 mm
Blade Length: 2.375" (6 cm)
Overall Length: 7.375" (18,7 cm)

5472-08 [8 mm Offset Cement Removal Chisel]

Osteotome Dimensions: 8 mm Wide x 12 mm Long
Blade Length: 2.375" (6 cm)
Overall Length: 7.375" (18,7 cm)


5474-06 [6 mm Notched Cement Removal Chisel]

Osteotome Width: 6 mm
Blade Length: 2.625" (6 cm)
Overall Length: 7.375" (18,7 cm)

5475-08 [8 mm Implant Remover]

Diameter: 8 mm
Blade Length: 2.625" (6 cm)
Overall Length: 7.375" (18,7 cm)

5470-CASE [Case Only]



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.