


Engl Intercondylar Notch Retractors

Designed by Gerard A. Engl, MD

Enhances minimally invasive exposure of the medial femoral condyle in unicondylar knee arthroplasty



PRODUCT NO'S:

3230-01 [Small]
 Blade Width at Teeth: 9 mm
 Depth from Bend: 2.25" (5.7 cm)
 Overall Length: 8.125" (20.6 cm)

3230-02 [Medium]
 Blade Width at Teeth: 10 mm
 Depth from Bend: 2.25" (5.7 cm)
 Overall Length: 8.125" (20.6 cm)

3230-03 [Large]
 Blade Width at Teeth: 12 mm
 Depth from Bend: 2.25" (5.7 cm)
 Overall Length: 8.125" (20.6 cm)



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
 info@innomed.net

912.236.0000 Phone
 912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
 info@innomed-europe.com

1.800.548.2362