

Evans Universal Carpal Tunnel Knife Guide

Designed by Peter J. Evans, MD, PhD

Designed to protect the median nerve while providing a choice of grooved tracks for commercially available retrograde knives (that do not provide this feature) or for tenotomy scissors

Allows for smooth advance of the blade or scissors to divide the transverse carpal ligament.
Designed for a mini-open, non-endoscopic approach.

PRODUCT NO:

1128

Overall Length: 8" (20,3 cm)

Blade Guide Widths: 2 mm and 5 mm



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362