


Extended Cup Positioner

Designed modification by James F. Kayanfar, MD of an original design by Thomas Eickmann, MD

Designed to help reposition an acetabular cup during total hip arthroplasty



Ultra hard titanium nitride coating helps to extend life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.

PRODUCT NO:

5475-10

Diameter: 8 mm

Overall Length: 12.75" (32,4 cm)

Handle Length: 4.75" (12,1 cm)

Shaft Length: 8" (20,3 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362