

Leg Stabilizer

Designed by Gregory Fanelli, MD

*Useful in arthroscopic knee surgery
to hold the leg in position*

Helps to open up the knee joint when pressure is applied to the lower leg. Pad and sterilizable table clamp included.

PRODUCT NO'S:

8840

Overall Length: 18.5" (47 cm)


Handle Length: 9.25" (23,5 cm)

Pad Diameter: 3" (7,6 cm)

Replacement Parts:

9120 [Table Clamp]

8840-P [Pad]



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362