

Femoral Tibial Spreader with Speed Lock

Designed by Nasim A. Rana, MD

Helps separate the femur and tibia during total knee replacement surgery



PRODUCT NO:

1855-SL [Speed Lock]
Overall Length: 10" (25,4cm)
Pads: 23x14mm
Opens to 50mm

**PROUDLY
MADE
IN THE
USA**

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362