


Femoral Tibial Coated Spreader Bar

Designed to separate the femur and tibia when implant components are in place

The end is coated to help protect from scratching component surfaces.

Product No: 1820
Overall Length: 13"
Coated Surface: 4"
Blade Width: 13mm



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

