


# Femoral Trial Extraction Hook

*Designed to help remove a trial femoral component during total knee arthroplasty*



Product No's:  
3635 [Hook with Slap Hammer]  
3635-01 [Hook Only]  
*Overall Length: 2.85"*  
*Hook End: 5mm X 5mm*  
3040 [Slap Hammer]  
*Overall Length: 6.75" - 10"*

**FREE TRIAL**  
ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**  
912 236 0000  
Fax 912 236 7766  
Internet [www.innomed.net](http://www.innomed.net)  
eMail [info@innomed.net](mailto:info@innomed.net)

**Innomed**  
ORTHOPEDIC INSTRUMENTS  
103 Estus Drive Savannah GA 31404

