

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2002 Innomed, Inc., All Rights Reserved.

Innomed

Instrument Division
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Femoral Canal Preparation Scrapers

Used to help remove excess bone in the inside of the femoral canal after broaching. Available in two blade widths.

Designed by
Gary W. Bradley, M.D.

Product No's: 5450-01 [9mm]
5450-02 [12mm]

*Used for femoral canal preparation
in total hip surgery*

