


Flexible Ball Nose Reamer

Designed by Stu Allen

Designed for safe and effective use in removing pedestal formation in the femoral and tibial canals

Recommended for use with a guide wire. Cannulated to allow guide wire use. Connects with a Zimmer Hall Quick-connect.



Quick-connect for use with a driver.

PRODUCT NO'S:

2628 [Flexible Ball Nose Reamer with Zimmer Hall Quick-connect]

Overall Length: 10" (25,4 cm)

Reamer Diameter: 7,5 mm

8248 [Fixed Driver with Zimmer Hall Quick-connect]

Overall Length: 5.75" (15,6 cm)

Handle Width: 4.625" (11,6 cm)



Fixed Driver
with Zimmer Hall
Quick-connect

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Scan to
Launch Our
Website



Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362