




# Gap Clamp for Cortical Button Distal Bicep Repair

Designed by Corey Trease, MD

*Designed to be used to help consistently set  
the gap for the radius cortex between the distal  
biceps stump and the cortical button*



**PRODUCT NO:**  
5262  
Overall Length: 8" (20.4 cm)  
Jaw Width Tapers from: 6 to 2.5 mm

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**