

Goldstein Flanged Femoral Elevator

Designed to elevate the femoral neck and prevent soft tissues from covering the medial femoral neck


Designed by
Wayne M. Goldstein, MD



Allows enhanced visualization of the neck osteotomy even in a small incision. The flange also helps hold the elevator securely under the neck so it will not slip out during the vibration with broaching.

Products No's:
7655-R [Right]
Overall Length: 20.25"
Handle Length: 12"
Blade Width: 38mm
Blade Flange: 10mm x 17mm

7655-L [Left]
Overall Length: 20.25"
Handle Length: 12"
Blade Width: 38mm
Blade Flange: 10mm x 17mm



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2009 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

