

# Goldstein Suture Cutter

Designed by Wayne M. Goldstein, MD

*Allows for the use of a replaceable standard #10 scalpel blade, ensuring a precise cut and reducing the risk of tearing a suture, suture material, or surrounding tissue*


The entire unit is autoclavable.

**PRODUCT NO:**

5159

Overall length: 8.25"

Shaft Length: 4.25"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

© 2011 Innomed, Inc.

INNOMED

1.800.548.2362