



# Gorski Hip Cup Extraction Hook

Designed by Jerrold Gorski, MD

*Helps in the removal of a hip cup*

Designed to quickly fit into a screw hole of a hip cup after the screws have been removed and the cup loosened. The slap hammer helps to remove the cup in the angle it was inserted.

## PRODUCT NO'S:

### Hook for 5.0mm Screw Holes

3665 [Hook with Standard (#3925) Slap Hammer]

3665-01 [Hook w/o Slap Hammer]

### Optional Items:

3935 [Extra Large Slap Hammer] Thread Gauge: 3/8"-16



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

# INNOMED