

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Gorski Hip Cup Extraction Hook

Designed by
Jerrold Gorski, MD

Designed to quickly fit into a screw hole of a hip cup after the screws have been removed and the cup loosened. The slap hammer helps to remove the cup in the angle it was inserted.

Product No's:
3660 [Hook for 6.5mm Screw Holes with Slap Hammer]
3660-01 [Hook for 6.5mm Screw Holes - No Slap Hammer]
3665 [Hook for 5.0mm Screw Holes with Slap Hammer]
3665-01 [Hook for 5.0mm Screw Holes - No Slap Hammer]
3935 [Slap Hammer]

*Helps in the
removal of a
hip cup*

