

Use punch to
remove tibial
bone plug, then...

...use tamp to insert
plug as autograft for the
femoral intramedullary
alignment hole



Goytia Osteotome Punch Tamp Assembly

Designed by Robin Goytia, MD

*Designed for removing a tibial bone plug to
use as autograft for the femoral intramedullary
alignment hole in total knee replacement*

PRODUCT NO'S:

5339-00 [Punch & Tamp Set]

Set Includes / Available Individually:

5339-01 [Osteotome Punch]

Overall Length: 7.75" (19,7 cm)



Outside Diameter: 16 mm

Inside Diameter: 13.7 mm

5339-02 [Tamp]

Overall Length: 7.75" (19,7 cm)

Diameter: 12.3 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362