





Grant TKA Anatomic Bone File Set

Designed by Richard E. Grant, MD

A bone rasp and plumb rod set designed for TKA tibial cut surface preparation




Patent Pending



Plumb rod fits into the handle of each bone rasp: 0°, 2° Left, and 2° Right.

PRODUCT NO'S:	
6906-00	[Set]
Set Includes/ Available Separately:	
6906-01	[Plumb Rod] Overall Length: 14" (35,6 cm)
6906-02	[0° (Flat) Rasp] Overall Length: 6.375" (16,2 cm) Rasp Platform Length: 1.7" (4,3 cm) Rasp Platform Width: 2.7" (6,9 cm)
6906-03	[2° Right Rasp] Overall Length: 6.375" (16,2 cm) Rasp Platform Length: 1.7" (4,3 cm) Rasp Platform Width: 2.7" (6,9 cm)
6906-04	[2° Left Rasp] Overall Length: 6.375" (16,2 cm) Rasp Platform Length: 1.7" (4,3 cm) Rasp Platform Width: 2.7" (6,9 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2021 Innomed, Inc.