



Prototype Shown



Prototype Shown



Prototype Shown

# Hand/Finger Positioner

Designed by Emad Aboujaoude, MS, MPAS, PA-C

*Designed to help provide surgical positioning during flouroscopy and fixation by isolating the operative digit while retracting the unaffected digits*

Radiolucent positioner can be steam or gas sterilized.

Uses include but not limited to:

- ▶ Intramedullary Metacarpal Screw
- ▶ Phalanges CRPP
- ▶ Digit Amputation
- ▶ Digit Mass Excision
- ▶ Finger Joint Fusion



**PRODUCT NO:**

1134

Overall Length: 5.75" (14,6 cm)  
Handle Width: 4.25" (10,8 cm)  
Blade Width: 1.6" (4 cm)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe**

Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**

# INNOMED