



Hannum Modified Angled Grasper

Designed by Scott Hannum, MD

Heavy duty large bone grasper designed to help trim acetabular osteophytes – angled to ergonomically fit around the rim via the direct anterior approach



PRODUCT NO:
1775-04
Overall Length: 8.5" (21,6 cm)
Jaw Width: 11 mm
Jaw Bite Internal: 9 mm x 21 mm



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71



1.800.548.2362