

Hannum Rongeur/Grasper

Designed by Scott Hannum, MD

Teeth in jaw firmly holds bone and tissue



Non-locking design can be easily gripped while allowing greater pressure to be applied. Available in three jaw sizes: short jaw for holding bone, medium jaw for smaller bones, and long jaw for tissue.

PRODUCT NO'S:

1775-01 [Short Jaw]
8mm Jaw Width
Overall Length: 9.25"

1775-02 [Medium Jaw]
5mm Jaw Width
Overall Length: 9.25"

1775-03 [Long Jaw]
3mm Jaw Width
Overall Length: 9.25"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362