

Havens Modified Kolbel Soft Tissue Retractor

Designed by Philip Havens, MD

Designed for retraction on deltoid split incisions on mini-open rotator cuff repairs



Jaws and arms are parallel with no gap when closed to allow easier insertion in tight spaces.

PRODUCT NO:

T1006-02

Overall Length: 7.5" (19,1 cm)

Opens To: 4.5" (11,4 cm)

Prong Depth: 18 mm



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362