

# Hibbs Retractors

*Designed for soft tissue retraction by either the toothed end or curved handle end*

Extra large used in large patients when more leverage and depth is needed.

Product No's:

6230 [Extra Large]  
Overall Length: 14.25"  
Handle Length: 13"  
Blade Depth: 6.5"  
Blade Width: 38mm

6235 [Medium]  
Overall Length: 10.75"  
Handle Length: 9.75"  
Blade Depth: 4.5"  
Blade Width: 25mm

6240 [Standard]  
Overall Length: 8.75"  
Handle Length: 8"  
Blade Depth: 3"  
Blade Width: 25mm



**FREE TRIAL**  
ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

**Phone** 800 548 2362  
912 236 0000  
**Fax** 912 236 7766  
**Internet** [www.innomed.net](http://www.innomed.net)  
**eMail** [info@innomed.net](mailto:info@innomed.net)

**Innomed**  
ORTHOPEDIC INSTRUMENTS  
103 Estus Drive Savannah GA 31404

