

Retractors for Hip Surgery – Double Bent

For general use in total, minimally invasive, and revision hip surgery

A Single Prong Double Bent Hohmann Acetabular Retractor

PRODUCT NO'S:

6210 [Standard Blade]
Overall Length: 10.5" (26,7 cm)
Blade + Tip Length: 2.5" (64 mm)
Blade Width: 15 mm

6212 [Longer Blade]
Overall Length: 11.25" (28,6 cm)
Blade + Tip Length: 3.5" (89 mm)
Blade Width: 15 mm



A Single Prong Double Bent Hohmann Acetabular Retractor

A Single & Double Prong Double Bent Hohmann Acetabular Retractor – Long

Extra Grip Tip design modification by Alfred A. Durham, MD

PRODUCT NO'S:

6210-02 [Single Prong Standard Blade]
Overall Length: 12.5" (31,8 cm)
Blade + Tip Length: 3" (76 mm)
Blade Width: 15 mm

6210-02L [Lighted Single Prong Standard Blade]
Overall Length: 12.5" (31,8 cm)
Blade + Tip Length: 3" (76 mm)
Blade Width: 15 mm

6211 [Single Prong Standard Blade with Extra Grip Tip]
Overall Length: 12.5" (31,8 cm)
Blade + Tip Length: 3" (76 mm)
Blade Width: 15 mm

6213 [Single Prong Longer Blade]
Overall Length: 15" (38,1 cm)
Blade + Tip Length: 5" (127 mm)
Blade Width: 15 mm

6220 [Double Prong Standard Blade]
Overall Length: 12.5" (31,8 cm)
Blade + Tip Length: 3" (76 mm)
Blade Width: 15 mm



Extra Grip Tip helps to stabilize the retractor around the bone

A Single & Double Prong Double Bent Hohmann Acetabular Retractor – Long

Lighted retractor can be attached to a fiber optic light cable with ACMI (female) connector.

A Single Prong Double Bent Hohmann Acetabular Retractor – Extra Long

PRODUCT NO'S:

6210-04 [Standard Blade]
Overall Length: 16.25" (41,3 cm)
Blade + Tip Length: 3" (76 mm)
Blade Width: 15 mm

6214 [Longer Blade]
Overall Length: 18" (45,7 cm)
Blade + Tip Length: 5" (127 mm)
Blade Width: 15 mm



A Single Prong Double Bent Hohmann Acetabular Retractor – Extra Long

www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS



ISO 13485:2016




103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362



Retractors for Hip Surgery – Curved

For general use in total, minimally invasive, and revision hip surgery



B Single Prong Broad Acetabular Retractor

Single Prong Broad
#6320



C Double Prong Broad Acetabular Retractor

Double Prong Broad
#6160



Standard
#6450



Standard with Short Tip
#6450-03



Extra Deep
#6450-01



Extra Deep with Short Tip
#6450-04



Straight Tip
#6450-02

C Single Prong Soft Tissue Retractors



Standard
#6570

D Single Prong Acetabular Retractors

Extra Deep
#6570-01

B Single Prong Broad Acetabular Retractor

PRODUCT NO:
6320
Overall Length: 12" (30,5 cm)
Blade Width: 40 mm

C Double Prong Broad Acetabular Retractor

PRODUCT NO:
6160
Overall Length: 12.5" (31,8 cm)
Blade Width: 40 mm

C Single Prong Soft Tissue Retractors

PRODUCT NO:
6450 [Standard]
Overall Length: 12.125" (30,8 cm)
Blade Width: 22.3 mm
6450-03 [Standard with Short Tip]
Overall Length: 11.75" (29,8 cm)
Blade Width: 22.3 mm
6450-01 [Extra Deep]
Overall Length: 13.75" (34,9 cm)
Blade Width: 22.3 mm
6450-04 [Extra Deep w/Short Tip]
Overall Length: 13" (33 cm)
Blade Width: 22.3 mm
6450-02 [Straight Tip]
Overall Length: 12.125" (30,8 cm)
Blade Width: 22.3 mm

D Single Prong Acetabular Retractors

PRODUCT NO:
6570 [Standard]
Overall Length: 12.125" (30,8 cm)
Blade Width: 22.3 mm
6570-01 [Extra Deep]
Overall Length: 13.75" (34,9 cm)
Blade Width: 22.3 mm