

Retractors for Hip Surgery – Double Bent

For general use in total, minimally invasive, and revision hip surgery



A Single Prong Double Bent Hohmann Acetabular Retractor

A Single Prong Double Bent Hohmann Acetabular Retractor

PRODUCT NO'S:
6210 [Standard Blade] Overall Length: 10.5" (26,7 cm) Blade + Tip Length: 2.5" (64 mm) Blade Width: 15 mm
6212 [Longer Blade] Overall Length: 11.25" (28,6 cm) Blade + Tip Length: 3.5" (89 mm) Blade Width: 15 mm

A Single & Double Prong Double Bent Hohmann Acetabular Retractor – Long

Extra Grip Tip design modification by Alfred A. Durham, MD

PRODUCT NO'S:
6210-02 [Single Prong Standard Blade] Overall Length: 12.5" (31,8 cm) Blade + Tip Length: 3" (76 mm) Blade Width: 15 mm
6210-02L-01 [Lighted Single Prong Standard Blade] Overall Length: 12.5" (31,8 cm) Blade + Tip Length: 3" (76 mm) Blade Width: 15 mm
6211 [Single Prong Standard Blade with Extra Grip Tip] Overall Length: 12.5" (31,8 cm) Blade + Tip Length: 3" (76 mm) Blade Width: 15 mm
6213 [Single Prong Longer Blade] Overall Length: 15" (38,1 cm) Blade + Tip Length: 5" (127 mm) Blade Width: 15 mm
6220 [Double Prong Standard Blade] Overall Length: 12.5" (31,8 cm) Blade + Tip Length: 3" (76 mm) Blade Width: 15 mm

Lighted retractor can be attached to a fiber optic light cable with ACMI (female) connector.

A Single Prong Double Bent Hohmann Acetabular Retractor – Extra Long

PRODUCT NO'S:
6210-04 [Standard Blade] Overall Length: 16.25" (41,3 cm) Blade + Tip Length: 3" (76 mm) Blade Width: 15 mm
6214 [Longer Blade] Overall Length: 18" (45,7 cm) Blade + Tip Length: 5" (127 mm) Blade Width: 15 mm



A Single & Double Prong Double Bent Hohmann Acetabular Retractor – Long

Extra Grip Tip helps to stabilize the retractor around the bone



A Single Prong Double Bent Hohmann Acetabular Retractor – Extra Long

Disposable LED Light Source

PRODUCT NO'S:
PACKAGE OF 1:
8010-01 [Disposable LED Light Source] Overall Length: 2.5" (6,4 cm) Diameter: 1" (2,54 cm)
PACKAGE OF 10:
8010-10 [Disposable LED Light Source]



Lighted retractor comes with one (1) Disposable LED Light Source (#8010-01)

ISO 13485:2016

Scan to Launch Our Website



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net


912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMEDI

1.800.548.2362





Retractors for Hip Surgery – Curved

For general use in total, minimally invasive, and revision hip surgery



B Single Prong Broad Acetabular Retractor

Single Prong Broad
#6320



C Double Prong Broad Acetabular Retractor

Double Prong Broad
#6160



Standard
#6450



Standard with Short Tip
#6450-03

New!



Extra Deep
#6450-01



Extra Deep with Short Tip
#6450-04

New!



Straight Tip
#6450-02

C Single Prong Soft Tissue Retractors



Standard
#6570

D Single Prong Acetabular Retractors

Extra Deep
#6570-01

B Single Prong Broad Acetabular Retractor

PRODUCT NO:
6320
Overall Length: 12" (30,5 cm)
Blade Width: 40 mm

C Double Prong Broad Acetabular Retractor

PRODUCT NO:
6160
Overall Length: 12.5" (31,8 cm)
Blade Width: 40 mm

C Single Prong Soft Tissue Retractors

PRODUCT NO:
6450 [Standard]
Overall Length: 12.125" (30,8 cm)
Blade Width: 22.3 mm
6450-03 [Standard with Short Tip]
Overall Length: 11.75" (29,8 cm)
Blade Width: 22.3 mm
6450-01 [Extra Deep]
Overall Length: 13.75" (34,9 cm)
Blade Width: 22.3 mm
6450-04 [Extra Deep w/Short Tip]
Overall Length: 13" (33 cm)
Blade Width: 22.3 mm
6450-02 [Straight Tip]
Overall Length: 12.125" (30,8 cm)
Blade Width: 22.3 mm

New!

New!

D Single Prong Acetabular Retractors

PRODUCT NO:
6570 [Standard]
Overall Length: 12.125" (30,8 cm)
Blade Width: 22.3 mm
6570-01 [Extra Deep]
Overall Length: 13.75" (34,9 cm)
Blade Width: 22.3 mm

ISO 13485:2016

Scan to
Launch Our
Website



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71



1.800.548.2362