


Horseshoe Shoulder Frame and Blade Assembly

Designed to enhance exposure during shoulder arthroplasty procedures



PRODUCT NO'S:

2030-00 [Set]
Includes (1) Frame, (1) Wide Blade, (1) Narrow Blade,
and (1) Narrow, Bent Blade


Also available individually:

2030-01 [Horseshoe Frame]
Overall Dimensions: 7" x 5" (17,8 cm x 12,7 cm)
Frame Width: .7" (15 mm)

2030-02 [Retractor Blade – Wide]
Blade Width: 22 mm
Overall Length: 11" (27,9 cm)

2030-03 [Retractor Blade – Narrow]
Blade Width: 14 mm
Overall Length: 11" (27,9 cm)

2030-04 [Retractor Blade – Narrow, Bent]
Blade Width: 14 mm
Overall Length: 10" (25,4 cm)
Handle Length: 4.5" (11,4 cm)



ISO 13485:2016

Scan to
Launch Our
Website



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362