

J.B. Redler Retractor

Designed by M.R. Redler, MD

*Uniquely balanced retractor for bone exposure
for a multitude of upper extremity procedures*

Double-angle design allows for ideal exposure with minimal effort to hold the retractor, while the assistant's hands are well out of the way of the exposure. The aperture in the base of the handle allows the retractor to be attached via a Penrose drain to the table for hands-free approach.

PRODUCT NO:

1645

Overall Length: 5"

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362