

Universal Table Adapter with Traction Device Assembly

A universal traction assembly with Jackson Traction Device for standard operating room tables used to assist with fracture fixation in the acetabulum, pelvis, and femur, and designed to attach to standard operating table side rails

PRODUCT NO'S:

0005-00 [Complete Assembly]
<i>Assembly Includes / Available Individually:</i>
0005-01 [Universal Table Adapter]
0005-02 [Vertical Rod] <i>Two included in set; One with this product number</i>
0005-03 [Horizontal Rod]
0007 [Jackson Traction Device] <i>This product number includes (1) #0008 Disposable Sterile Kit</i>
Optional Items:
0008-CASE [Case of Sterile Kits] Pkg of 10
2595 [Table Clamp]
Replacement Item:
4150-PS [Post Screw]



All components can be ordered separately.

PRODUCT NO'S:

0005-10 [Assembly without #0007 Jackson Traction Device]
<i>Assembly Includes / Available Individually:</i>
0005-01 [Universal Table Adapter]
0005-02 [Vertical Rod] <i>Two included in set; One with this product number</i>
0005-03 [Horizontal Rod]
Optional Item:
2595 [Table Clamp]
Replacement Item:
4150-PS [Post Screw]

