



# Joint, Calcaneal, Small Bone Compressor/Distractor

*Selection lever switches the mechanism from compression to distraction*

Simply squeeze the handle one time after direction selection to engage the mechanism.  
Two hole sizes for pin size selection. Calibrated ratchet handle.



#### PRODUCT NO'S:

4865-LS [Standard] Overall Length: 8.5" Holes For: .062" & .094" (1.6 & 2.5 mm) K-wire Pins	 <p>Thumbscrews help prevent the unit from sliding on the pins.</p>
4865-LS-TS [With Thumbscrews] Overall Length: 8.5" Holes For: .062" & .094" (1.6 & 2.5 mm) K-wire Pins	

ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**