

Kawell Short Army Navy Retractor

Designed by Ron Kane, DPM

A short handled Army Navy retractor, especially useful with a gastrocnemius recession



PRODUCT NO:

1148

Overall Length: 4.75" (12.1 cm)

Large End Blade Length: 1.75" (4.4 cm)

Large End Blade Width: .625" (1.6 cm)

Small End Blade Length: 1" (2.5 cm)

Small End Blade Width: .3125" (0.8 cm)



USA MADE

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362