


Kenerly Double Parallel Femoral Neck Disc Grasper

Designed by J. Lex Kenerly, III, MD

Designed to remove the central disc of a double, parallel cut femoral neck osteotomy when performing THA





PRODUCT NO:

1720-01

Overall Length: 10" (25,4 cm)

Jaw Dimensions: 1.44" x .72" (36,6 cc x 18,3 mm)

Lower Jaw Thickness: 1 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362