


Kolbel Self-Retaining Glenoid Retractors



**Radiolucent
Blades
Now Available!**

The OrthoLucent™ carbon fiber PEEK blade is strong, lightweight, completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.

Kolbel Self-Retaining Retractor Blades



PRODUCT NO'S:	
Wide Blades	
T1018	[36 x 36 mm]
T1019	[36 x 53 mm]
T1020	[36 x 68 mm]
T1021	[36 x 85 mm]
Radiolucent Blade	
T1019-R*	[36 x 53 mm]

Narrow Blades	
T1022	[20 x 36 mm]
T1023	[20 x 53 mm]
T1024	[20 x 68 mm]
T1025	[20 x 85 mm]



Kolbel Self-Retaining Glenoid Retractor with Hinge and Ergonomic Handle

Designed with longer articulating arms—helpful for use with larger patients

Two pairs of snap-in, freely pivoting blades included.

PRODUCT NO'S:	
T1016-01	[Set]
Set Includes:	
T1016-01-2F	[Retractor] Overall Length: 10.75" (27,3 cm) Length-to-hinge: 7.75" (19,7 cm) Arm Length: 3 (7,6 cm)
T1018-P	[Blades—Pair] 36 mm X 36 mm
T1019-P	[Blades—Pair] 36 mm X 53 mm




Double Ergonomic Handle
For added control and comfort

Modified Kolbel Self-Retaining Glenoid Retractor with Hinge

Two pairs of snap-in, freely pivoting blades included.

PRODUCT NO'S:	
T1014-01	[Set – Standard Handle]
T1014-01-2F	[Set – Ergonomic Handle]
Set Includes:	
T1015-01	[Retractor – Standard Handle] Overall Length: 8.25" (21 cm) Length-to-hinge: 6" (15,2 cm) Arm Length: 2.25 (5,7 cm)
– OR –	
T1015-01-2F	[Retractor – Ergonomic Handle] Overall Length: 9.25" (23,5 cm) Length-to-hinge: 7" (17,8 cm) Arm Length: 2.25 (5,7 cm)
T1018-P	[Blades—Pair] 36 mm X 36 mm
T1019-P	[Blades—Pair] 36 mm X 53 mm




Standard Handle

Ergonomic Handle
For added control and comfort

Kolbel Self-Retaining Glenoid Retractor

Two pairs of snap-in, freely pivoting blades included.

PRODUCT NO'S:	
T1014	[Set – Standard Handle]
T1014-2F	[Set – Ergonomic Handle]
Set Includes:	
T1015	[Retractor – Standard Handle] Overall Length: 8.25" (21 cm)
– OR –	
T1015-2F	[Retractor – Ergonomic Handle] Overall Length: 9.25" (23,5 cm)
T1018-P	[Blades—Pair] 36 mm X 36 mm
T1019-P	[Blades—Pair] 36 mm X 53 mm



Standard Handle

Ergonomic Handle
For added control and comfort



– See Other Side for More Kolbel Retractors –

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362



Kolbel Self-Retaining Glenoid Retractors



**Radiolucent
Blades
Now Available!**

The OrthoLucent™ carbon fiber PEEK blade is strong, lightweight, completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.

Kolbel Self-Retaining Retractor Blades



PRODUCT NO'S:	
Wide Blades	
T1018	[36 x 36 mm]
T1019	[36 x 53 mm]
T1020	[36 x 68 mm]
T1021	[36 x 85 mm]
Radiolucent Blade	
T1019-R*	[36 x 53 mm]
Narrow Blades	
T1022	[20 x 36 mm]
T1023	[20 x 53 mm]
T1024	[20 x 68 mm]
T1025	[20 x 85 mm]




* MADE EXCLUSIVELY FOR INNOMED IN SWITZERLAND

Kolbel Self-Retaining Glenoid Retractor with Center Blade

Center blade can be reversed for shallow or deep retraction

Two pairs of snap-in, freely pivoting blades included.

PRODUCT NO'S:	
T1050	[Set – Standard Handle]
T1050-2F	[Set – Ergonomic Handle]
Set Includes:	
T1050-01	[Retractor – Standard Handle] Overall Length: 8" (20,3 cm)
– OR –	
T1050-01-2F	[Retractor – Ergonomic Handle] Overall Length: 9" (22,9 cm)
T1050-02	[Center Blade] Length-to-bend: 6.25" (15,9 cm) Depth: 2.5" (6,4 cm)
T1018-P	[Blades–Pair] 36 mm X 36 mm
T1019-P	[Blades–Pair] 36 mm X 53 mm



Standard Handle

Ergonomic Handle
For added control and comfort

Kolbel Self-Retaining Retractor

Two pairs of snap-in, freely pivoting blades included.

PRODUCT NO'S:	
T1016	[Set]
Set Includes:	
T1017	[Retractor] Overall Length: 8.25" (21 cm) Arm Length: 6.125" (15,6 cm) Arm Length-to-hinge: 3" (7,6 cm)
T1018-P	[Blades–Pair] 36 mm X 36 mm
T1019-P	[Blades–Pair] 36 mm X 53 mm



– See Other Side for More Kolbel Retractors –

ISO 13485:2016

Scan to
Launch Our
Website



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362