

Kudrna Cup Channel Chisel

Designed by James C. Kudrna, MD

Designed to help break the bone-prosthetic interface of well-fixed non-cemented acetabular components being revised

Ultra hard titanium nitride coating helps to extend chisel life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.

PRODUCT NO:

3686

Overall Length: 12" (30,5cm)

Handle Length: 3.5" (8,9cm)

Blade Width: 30mm

Blade Depth: 15mm

**PROUDLY
MADE
IN THE
USA**



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED The Innomed logo features the word 'INNOMED' in a bold, sans-serif font, followed by a stylized lowercase 'i' that incorporates a human figure silhouette. To the right of the logo is a small American flag.

1.800.548.2362