


Calcaneal Lateral Column Spreader

Designed by Keith L. Wapner, MD

Used for lateral column lengthening of the calcaneus

PRODUCT NO:

1725
Pads: 14mm x 12mm
Arms Open to: 45mm
Overall Length: 4.25"



FREE TRIAL ON MOST INSTRUMENTS

© 2011 Innomed, Inc.

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

1.800.548.2362