



# Lawton Distal Radius Mini Frame & Blade Set


Designed by Jeffrey Lawton, MD

*Designed for self-retaining exposure for distal radius and other small bone fractures*



Short Blade Depth

Small & Large Blade Depth



Short Blade

Small Blade

Optional Large Blade

Set includes: (1) Frame, (2) Short Blades, (2) Small Blades.  
Optional Large Blade available separately.

PRODUCT NO:
<b>1578-00</b> [Set]
<b>Set Includes / Available Individually:</b>
<b>1578-01</b> [Mini Frame] <i>Dimensions: 3" x 2.5" (7,6 x 6,4 cm)</i>
<b>1578-02</b> [Mini Short Blade] (2) included in set, (1) with this product number <i>Overall Length: 2.5" (6,4 cm)</i> <i>Blade Width: .625" (16 mm)</i> <i>Blade Depth: .875" (22 mm)</i>
<b>1578-03</b> [Mini Small Blade] (2) included in set, (1) with this product number <i>Overall Length: 2.625" (6,7 cm)</i> <i>Blade Width: .625" (16 mm)</i> <i>Blade Depth: 1.125" (29 mm)</i>
<b>Optional Blade / Not Included In Set:</b>
<b>1578-04</b> [Mini Large Blade] <i>Overall Length: 2.5" (6,4 cm)</i> <i>Blade Width: .935" (24 mm)</i> <i>Blade Depth: 1.125" (29 mm)</i>



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362