



Lee MIS Bunion Parallel Guide

Designed by Wonyong Lee, MD

Designed to facilitate MIS bunion correction, focusing on achieving the ideal parallel trajectory of two guide pins for 1st metatarsal screw fixation, with a beveled guide hole to help prevent pin slippage off the bone



PRODUCT NO:
3023
 Overall Length: 6.375 (17.1 cm)
 Top Tube Guide Length: 1.375" (3.5 cm)
 Top Tube Hole Size: For 1,6 mm Pins
 Lower Tube Guide Length: 2.375" (6 cm)
 Lower Tube Hole Size: For 1,4 mm Pins
 Distance Between Holes: 6 mm



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED  

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
 info@innomed.net

912.236.0000 Phone
 912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
 info@innomed-europe.com

1.800.548.2362

© 2025 Innomed, Inc.